

V-32, Multi-Vortex



DESCRIPTION

Multi-Vortex **V-32** is intended for intensive stirring of bacterial and yeast cell, washing from the culture medium and extraction of metabolites and enzymes from cells and cell cultures.

It is different from V-1 by the possibility of mixing up to 32 tubes simultaneously.

Vortex is applicable for:

- Performing various DNA operations — deproteinisation of DNA/protein complexes;
- Purification of low-molecular DNA/RNA fragments in PCR-diagnostic.

Multi-Vortex has two operation modes:

1. Continuous operation;
2. Impulse operation.

V-32 is supplied with a 32-socket universal platform for Eppendorf type tubes up to 1.5 ml (1.5/0.5/0.2 ml — 16/8/8 sockets) **PV-32** and **PL-1** platform for mixing single tube up to 50 ml.

An optional 6-socket platform **PV-6/10** for 10 ml tubes (max. tube diameter 15 mm) or **PV-48** for 48 - 0.2 ml microtubes can be supplied on request.



SPECIFICATION

Speed control range	500-3000 rpm
Acceleration time	3 s
Orbit	2 mm
Continuous / impulse operation	+
Maximum load	70 g
Maximum continuous operation time	8 h
Dimensions	120x180x100 mm
Weight	1.5 kg
Input current/power consumption	12 V, 320 mA / 3.8 W
External power supply	Input AC 100-240 V; 50/60 Hz; Output DC 12 V



CAT. NR.

	V-32 incl. two platforms PV-32, PL-1 & external power supply unit
BS-010207-AAG	230VAC 50/60Hz Euro plug
BS-010207-AAK	230VAC 50/60Hz UK plug
BS-010207-AAK	230VAC 50/60Hz AU plug
BS-010207-AAK	100VAC 50/60Hz US plug
BS-010207-AAK	120VAC 60Hz US plug
-----	Optional accessories
BS-010207-BK	PV-6/10
BS-010207-GK	PV-48

BS-010207-EK	IQ OQ document
BS-010207-FK	PQ document

V-32, Multi-Vortex



-----	Replacement parts
BS-010207-CK	PV-32
BS-010207-AK	PL-1