
PSU-20i, Multi-functional Orbital Shaker





DESCRIPTION

Shaker **PSU-20i** provides three motion types: orbital, reciprocal and vibrating, which can be performed separately, pairwise and sequentially in repeated cycles. Shaker is designed for applications both in small specialized laboratories and in large multidisciplinary laboratories. **PSU-20i** is an ideal instrument for laboratories conducting research in biopharmaceutics and biomedicine.



Shaker PSU-20i is noiseless and reliable in operation, incorporates a direct drive system and brushless motor with a guaranteed service life up to 35,000 working hours. The use of direct drive and brushless motor allows for continuous mixing up to 7 days and ensures reliable operation for more than 2 years.

A choice of nine (9) different interchangeable platforms provides possibility of performing various procedures and techniques. Special attention should be paid to a multilevel platform, which allows accommodation of a large number of various microplates, Petri dishes, cultural bags and other low containers.

Shaker can be used in cold rooms or incubators, operating at ambient temperature range +4°C to +40°C.



Metion types	Description	Speed range	Turning angle	Motion timer *	General operation time (1 min step)
0 ⇔ orbital	Simple arbital reution with an op- tion of shifting direction (clockwise) anticlockwise) after set time	20-250 rpm	-	0-250 sec	1 min-96 hrs
0 ⇔ Reciprocal	Orbital motion with shifting direction of rotation	20-250 rpm	(10° increment)	0-250 sec	1 min-96 hrs
0 ∰ Vibrating	High speed, low amplitude motion	-	0°-5" (1" increment)	0-5 sec	1 min-96 hrs

^{*} for switching to the next motion in the cycle

SPECIFICATION

Speed control range	20-250* RPM (increment 5 rpm) * max. speed depends on the load and vessels' shape
Digital speed control	+
Maximum continuous operation time	168 hrs
Orbit	20 mm
Digital time setting	1 min-96 hrs / non-stop (increment 1 min)
Maximum load	8 kg
Overall dimensions (W \times D \times H)	410x410x130 mm
Weight	11.7 kg
Input current/power consumption	12 V, 3.2 A / 40 W
External power supply	Input AC 100–240 V; 50/60 Hz; Output DC 12 V

CAT. NR.

	PSU-20i with external power supply unit w-out platform
BS-010145-ACI	230VAC 50/60Hz Euro plug
BS-010145-ACQ	230VAC 50/60Hz UK plug
BS-010145-AC4	230VAC 50/60Hz AU plug
BS-010145-ACJ	100VAC 50/60Hz US plug
BS-010145-ACJ	120VAC 60Hz US plug
•	120 The Goriz Go plag
	Optional accessories
•	, ,
	Optional accessories
BS-010135-DK	Optional accessories Platform P-6/1000
BS-010135-DK BS-010135-AK	Optional accessories Platform P-6/1000 Platform P-9/500

PSU-20i, Multi-functional Orbital Shaker





BS-010135-JK	Platform UP-168
BS-010126-HK	FC-100, clamp 100ml
BS-010126-JK	FC-250, clamp 250ml
BS-010126-LK	FC-500, clamp 500ml
BS-010126-IK	FC-1000, clamp 1000ml
BS-010126-NK	FC-2000, clamp 2000ml
BS-010145-AK	Platform UP-330
BS-010145-BK	HB-330 additional holding bar for UP- 330
BS-010126-BK	Platform PP-20
BS-010126-BK BS-010126-CK	Platform PP-20 Platform PP-20-2
BS-010126-CK	Platform PP-20-2
BS-010126-CK BS-010126-DK	Platform PP-20-2 Platform PP-20-3
BS-010126-CK BS-010126-DK	Platform PP-20-2 Platform PP-20-3
BS-010126-CK BS-010126-DK BS-010126-EK	Platform PP-20-2 Platform PP-20-3 Platform PP-20-4

ACCESSORIES



P-6/1000

Platform with 6 clamps for 1000 ml flasks (360x400 mm)



P-9/500

Platform with 9 clamps for 500 ml flasks (360x400 mm)

PSU-20i, Multi-functional Orbital Shaker







P-16/250

Platform with 16 clamps for 250 ml flasks (360x400 mm)



P-30/100

Platform with 30 clamps for 100 ml flasks (360x400 mm)



UP-168

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately. (360x400 mm)



FC-100, FC-250, FC-500, FC-1000, FC-2000

Clamp for 100 ml - Ø65 mm, 250 ml - Ø85 mm, 500 ml - Ø105 mm, 1000 ml - Ø130 mm, 2000 ml - Ø166 mm flask (for UP-168)



UP-330

Universal platform with adjustable bars can accommodate laboratory glassware of different shapes $(345 \times 430 \times 105 \text{ mm})$

Biosan SIA., Ratsupites 7, build. 2, Riga, LV-1067, Latvia, Phone: +371 674 261 37, Fax: +371 674 281 01, E-mail: marketing@biosan.lv, http://www.biosan.ly

PSU-20i, Multi-functional Orbital Shaker







HB-330 Additional holding bar for UP-330



PP-20 --- PP-20-4

Flat platform with non-slip rubber mat can accommodate various low profile containers $(380 \times 480 \times 170/340/510 \text{ mm})$