

HiPo MPP-96, Microplate Photometer



DESCRIPTION

Microplate Photometer HiPo is a compact tabletop device for measuring the results of ELISA and microbiological studies in 96-well microplates. Photometer is controlled and outputs data via computer. The device is supplied with specialized software QuantAssay.

Features of QuantAssay software:

- ELISA assays of any complexity can be carried out via robust assay editor with help of Assay Wizard
- Quantitative assay includes up to 20 standards
- Avidity/Affinity assays
- Multiplex assays with up to 7 assays on one plate
- Qualitative assay includes up to 11 controls
- BestFit function for selecting the best calibration curve
- User friendly interface: get your results in 3 clicks
- Save, load and export results
- Creates visual reports

Accuracy (405, 450, 492, 620 nm)

0.000 - 2.000 OD $\leq (0.5 \% \pm 0.010 \text{ OD})$ typical

2.000 - 3.000 OD $\leq (1 \% \pm 0.010 \text{ OD})$ typical

Precision / Reproducibility (405, 450, 492, 620 nm)

0.000 - 2.000 OD $\leq (0.5 \% \pm 0.005 \text{ OD})$

2.000 - 3.000 OD $\leq (1.0 \% \pm 0.005 \text{ OD})$

* - up to 4 custom filters can be fitted on demand



SPECIFICATION

Detection mode	Absorbance
Light source	LED, self-calibrating
Photodetector	8 silicon photodiodes
Plate type	96-well microplates (including strip-well microplates)
Reading speed	5-8 s per wavelength
Measurement modes	Endpoint, kinetic and multi-label measurements
Measurement channels	8
Reference channel	1
Measurement range	0 - 4.3 OD
Resolution	0.0001 OD
Wavelength range	400 - 700 nm
Wavelength selection	up to 8* filters on wheel standard filters 405, 450, 492 and 620 nm

CAT. NR.

	HiPo MPP-96 incl. external power supply unit
BS-050108-A02	230VAC 50/60Hz Euro plug
BS-050108-A03	230VAC 50/60Hz UK plug
BS-050108-A05	100VAC 50/60Hz US plug
BS-050108-A05	120VAC 60Hz US plug
-----	Optional accessories
BS-050108-AK	<u>Verification instrument</u>
-----	Light Filters
400	400 nm
455	455 nm
458	458 nm
460	460 nm

HiPo MPP-96, Microplate Photometer



Shaking speed control range	4 amplitudes, 4 speeds
Software	QuantAssay
PC system requirements:	Intel/AMD Processor, 1 GB RAM, Windows Vista/7/8, USB
Weight	4.6 kg
Dimensions (WxDxH)	140 x 300 x 130 mm
External power supply	Input AC 100-240 V 50/60 Hz, Output DC 12 V

470	470 nm
480	480 nm
486	486 nm
488	488 nm
500	500 nm
508	508 nm
510	510 nm
515	515 nm
520	520 nm
532	532 nm
535	535 nm
540	540 nm
546	546 nm
550	550 nm
560	560 nm
568	568 nm
580	580 nm
589	589 nm
594	594 nm
600	600 nm
610	610 nm
632	632 nm
636	636 nm
640	640 nm
647	647 nm
650	650 nm
656	656 nm
660	660 nm
671	671 nm
676	676 nm
680	680 nm
685	685 nm
690	690 nm
694	694 nm